

NuBORNE™

Infant Warmer 500



Infant Radiant Warmer

- Control modes
 - Servo mode
 - Manual mode
 - Kangaroo mode
- LCD color display
- APGAR timer
- Vertical height adjustment
- Made in the USA

Warmology™ /wôrm-ol-o-gy/

The science of warm for the miracle of life

Features

NuBORNE™ Infant Warmer 500



Accessory Rail

Ideal for pole-mounted accessories including IV pole.



X-Ray Tray

Allows clinicians to take x-rays without disturbing the patient.



Instrument Shelf

Perfect for support equipment such as ventilators and cardiac monitors.



Patient Platform

Easy-open and removable sidewalls for access and cleaning.

Continuous Mattress Tilt: $\pm 12^\circ$
Mattress Size: 45cm x 65cm x 2.5cm



Storage Drawer

Fully removable and easy to clean storage drawer.



Vertical Height Adjustment

Quiet and smooth operation for an adjustable mattress height of 90cm to 110cm (floor to bed).

Order Information

Part #	Description
110-3395	NuBorne 500 Infant Warmer/220 VAC
108-4855	Temperature Probe (reusable)
110-3787	Temperature Probe (disposable) (10pk)
110-3833	IV Pole
110-3834	Monitor Stand
110-3836	Oxygen Cylinder Holder (short)
100-3835	Oxygen Cylinder Holder (long)



NICU

Europe

International Biomedical AG

Glütschbachstrasse 2, CH-3661 Uetendorf, Switzerland
T. +41.33.345.66.00 F. +41.33.345.66.08
welcome@int-bio.ch


Worldwide

International Biomedical

8206 Cross Park Drive, Austin, Texas USA 78754
T. 512.873.0033 F. 512.873.9090
sales@int-bio.com

int-bio.com

International Biomedical reserves the right to make changes in specifications and features shown herein or discontinue the product described at any time without notice or obligation. Contact your representative for the most current information.



Weight - 105kg
(Approx. Without Accessories)